




Product Datasheet

IDH1 antibody (orb37973)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	IDH1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Concentration	Concentration: Lot specific, refer to datasheet supplied with delivered product.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46659
Uniprot ID	O75874
NCBI	NP_001269316.1 , NP_001269315.1 , NP_005887.2
Dilution Range	IHC-P: 1:50-100, WB: 1:1000, IF/ICC: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of IDH1 antibody i...



Immunohistochemical analysis of formali...