

Product Datasheet

Phospho-GRK 2 (S685) Polyclonal Antibody (orb379730)

Description	Phospho-GRK 2 (S685) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P25098
Expiration Date	12 months from date of receipt.