



Product Datasheet

KHK Isoform C Antibody (orb379732)

Description	KHK Isoform C Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Concentration	1.0mg/mL
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P50053
Dilution Range	Western blotting: 1:500-1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of tumors d...



Western blot analysis of extracts from e...