


Product Datasheet

ST6GALNAC2 antibody (orb38044)

Description	ST6GALNAC2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 2 protein (101-374AA)
Target	ST6GALNAC2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UJ37
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 whole cell lysate.