

Product Datasheet

Muscle actin Mouse Monoclonal Antibody(2H3) (orb380515)

Description	Muscle actin Mouse Monoclonal Antibody(2H3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P68133
Expiration Date	12 months from date of receipt.