

Product Datasheet

Beclin Mouse Monoclonal Antibody(5C5) (orb380539)

Description	Beclin Mouse Monoclonal Antibody(5C5)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q14457
Expiration Date	12 months from date of receipt.