

Product Datasheet

MMP2 Mouse Monoclonal Antibody(1H1) (orb380542)

Description	MMP2 Mouse Monoclonal Antibody(1H1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P08253
Expiration Date	12 months from date of receipt.