

Product Datasheet

CHOP Mouse Monoclonal Antibody(2B1) (orb380543)

Description	CHOP Mouse Monoclonal Antibody(2B1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P35638
Expiration Date	12 months from date of receipt.