

Product Datasheet

ATG7 Mouse Monoclonal Antibody(3D6) (orb380613)

Description	ATG7 Mouse Monoclonal Antibody(3D6)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	095352
Expiration Date	12 months from date of receipt.