

Product Datasheet

Human Osteoprotegerin Fc Chimera protein (orb381018)

Description	Recombinant human Osteoprotegerin Fc Chimera protein
Reactivity	Human
Conjugation	Unconjugated
Target	TNFRSF11B
Preservatives	Lyophilized: 0.2 µm filtered concentrated solution in 20 mM PB, pH 6.0, 150 mM NaCl, 0.02 % Tween- 80
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 95% by SDS-PAGE and HPLC analyses.
Source	Yeast
MW	17.2 kDa
Expiration Date	6 months from date of receipt.