

Product Datasheet

CCDC6 Antibody (orb381052)

Description

CCDC6 Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, WB

Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of human CCDC6 (156-198aa KAELEQHLEQEQEFQVNLMMKKIKKLENDTISKQLTLEQLRR E), identical to the related mouse sequence.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

53291 MW

Uniprot ID


[Q16204](#)

Dilution Range


Western blot, 0.1-0.5µg/ml, Human, Rat
Immunocytochemistry/Immunofluorescence, 5µg/ml, Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Flow Cytometry analysis of CACO-2 cells ...



WB analysis of CCDC6 using anti-CCDC6 an...



IF analysis of CCDC6 using anti-CCDC6 an...