


Product Datasheet

Glucose 6 phosphate isomerase/GPI Antibody (orb381070)

Description

Glucose 6 phosphate isomerase/GPI Antibody



Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

WB analysis
of GPI using
anti-GPI
antibo...

Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of human GPI (5-39aa TRDPQFQKLQQWYREHRSELNLRRLFDANKDRFNH), different from the related mouse and rat sequences by sixteen amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

63147 MW

Uniprot ID

[P06744](#)

Dilution Range

Western blot, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.