


Product Datasheet

Human FGF1/Fgf Acidic ELISA Kit (orb381110)

Description	Human FGF1/Fgf Acidic ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A2-D155
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	17460 MW
Uniprot ID	P05230
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 200-400 pg/ml.
Expiration Date	Please enquire.



Human
FGF1
PicoKine
ELISA Kit
standard c...