


Product Datasheet

Mouse LIF ELISA Kit (orb381111)

Description	Mouse LIF ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S24-F203
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	22287 MW
Uniprot ID	P09056
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 7.8 pg/ml.
Expiration Date	Please enquire.



Mouse LIF
PicoKine
ELISA Kit
standard
CU...