


Product Datasheet

Mouse SPINT1/HAI-1 ELISA Kit (orb381115)

Description	Mouse SPINT1/HAI-1 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E30-L507
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	56591 MW
Uniprot ID	Q9R097
Sensitivity	12 pg/ml
Expiration Date	Please enquire.



Mouse
SPINT1/HAI-
1 PicoKine
ELISA Kit
st...