


Product Datasheet

Human HAI-2/SPINT2 ELISA Kit (orb381116)

Description	Human HAI-2/SPINT2 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M1-K197
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and cell lysates.
MW	28228 MW
Uniprot ID	043291
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human HAI-2/SPINT2 PicoKine ELISA Kit st...