


Product Datasheet

Mouse HAI-2/SPINT2 ELISA Kit (orb381117)

Description	Mouse HAI-2/SPINT2 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: R30-K140
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and cell lysates.
MW	27914 MW
Uniprot ID	Q9WU03
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 100 pg/ml.
Expiration Date	Please enquire.



Mouse HAI-2/SPINT2 PicoKine ELISA Kit st...