


Product Datasheet

Human TrkB ELISA Kit (orb381119)

Description	Human TrkB ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C32-H430
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and cell lysates.
MW	91999 MW
Uniprot ID	Q16620
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is 800-1000pg/ml
Expiration Date	Please enquire.



Human TrkB
PicoKine
ELISA Kit
standard c...