


Product Datasheet

Rat Hemojuvelin/RGM-C/HJV ELISA Kit (orb381123)

Description	Rat Hemojuvelin/RGM-C/HJV ELISA Kit
Reactivity	Rat
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q36-D396
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	45221 MW
Uniprot ID	Q8N7M5
Sensitivity	12 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 600-1000 pg/ml.
Expiration Date	Please enquire.



Rat
Hemojuvelin/RGM-
C PicoKine ELISA
Kit...