


Product Datasheet

Human FABP2/I-FABP ELISA Kit (orb381125)

Description	Human FABP2/I-FABP ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A2-D132
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	15207 MW
Uniprot ID	P12104
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 1-10 ng/ml.
Expiration Date	Please enquire.



Human
FABP2/I-
FABP
PicoKine
ELISA Kit
st...