


Product Datasheet

Human CD28 ELISA Kit (orb381128)

Description	Human CD28 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: N19-P152
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	25066 MW
Uniprot ID	P10747
Sensitivity	12 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 7.8 pg/ml.
Expiration Date	Please enquire.



Human CD28
PicoKine
ELISA Kit
standard c...