


Product Datasheet

Mouse TPOR/MPL/Thrombopoietin R ELISA Kit (orb381129)

Description	Mouse TPOR/MPL/Thrombopoietin R ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q26-W482
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	69788 MW
Uniprot ID	Q08351
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 156 pg/ml.
Expiration Date	Please enquire.



Mouse
TPOR/MPL
PicoKine
ELISA Kit
standa...