


Product Datasheet

Mouse EPCR ELISA Kit (orb381130)

Description	Mouse EPCR ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L18-S214
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27189 MW
Uniprot ID	Q64695
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:20, concentration in serum and plasma is around 20-30 ng/ml.
Expiration Date	Please enquire.



Mouse EPCR
PicoKine
ELISA Kit
standard c...