


Product Datasheet

Mouse FCRN/FCGRT ELISA Kit (orb381131)

Description	Mouse FCRN/FCGRT ELISA Kit
Reactivity	Mouse
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S22-S297
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and cell lysates.
MW	40093 MW
Uniprot ID	Q61559
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is 100-1000p/ml
Expiration Date	Please enquire.



Mouse
FCRN/FCGRT
PicoKine
ELISA Kit
stan...