


## Product Datasheet

### Mouse FABP4/A Fabp ELISA Kit (orb381132)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse FABP4/A Fabp ELISA Kit                                                                                                                                                      |
| <b>Reactivity</b>          | Mouse                                                                                                                                                                             |
| <b>Range</b>               | 156 pg/ml - 10,000 pg/ml                                                                                                                                                          |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                             |
| <b>Immunogen</b>           | Expression system for standard: E.coli; Immunogen sequence: M1-A132                                                                                                               |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Sample Types</b>        | cell culture supernatants, serum, plasma (heparin, EDTA) and cell lysates.                                                                                                        |
| <b>MW</b>                  | 14650 MW                                                                                                                                                                          |
| <b>Uniprot ID</b>          | <b>P04117</b>                                                                                                                                                                     |
| <b>Sensitivity</b>         | 15 pg/ml                                                                                                                                                                          |
| <b>Dilution Range</b>      | Dilution ratio of 1:50, concentration in serum and plasma is 20-100 ng/ml.                                                                                                        |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                   |



Mouse  
FABP4  
PicoKine  
ELISA Kit  
standard ...