


Product Datasheet

Human FABP4/A Fabp ELISA Kit (orb381133)

Description	Human FABP4/A Fabp ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-A132
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and cell lysates.
MW	14719 MW
Uniprot ID	P15090
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 5-20 ng/ml.
Expiration Date	Please enquire.



Human
FABP4
PicoKine
ELISA Kit
standard ...