


Product Datasheet

Rat MMP-9 ELISA Kit (orb381134)

Description	Rat MMP-9 ELISA Kit
Reactivity	Rat
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L35-P708
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	78611 MW
Uniprot ID	P50282
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is around 20 ng/ml.
Expiration Date	Please enquire.



Rat MMP-9
PicoKine
ELISA Kit
standard
CU...