


Product Datasheet

Human GHR/Growth Hormone R ELISA Kit (orb381136)

Description	Human GHR/Growth Hormone R ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A27-Y264
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	71500 MW
Uniprot ID	P10912
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 4 ng/ml.
Expiration Date	Please enquire.



Human
GHR
PicoKine
ELISA Kit
standard
cu...