


Product Datasheet

Human Fibulin-3/EFEMP1 ELISA Kit (orb381137)

Description	Human Fibulin-3/EFEMP1 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q18-F493
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	54641 MW
Uniprot ID	Q12805
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1000, concentration in serum is 1500-2800ng/ml
Expiration Date	Please enquire.



Human
Fibulin-
3/EFEMP1
PicoKine
ELISA Ki...