


Product Datasheet

Human HE4 ELISA Kit (orb381138)

Description	Human HE4 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: E31-F124
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.
MW	12993 MW
Uniprot ID	Q14508
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is 31.4-41.1 ng/ml.
Expiration Date	Please enquire.



Human HE4
PicoKine
ELISA Kit
standard
CU...