


Product Datasheet

Mouse Marapsin/Pancrestin ELISA Kit (orb381142)

Description	Mouse Marapsin/Pancrestin ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M38-V328
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	35789 MW
Uniprot ID	Q8BJR6
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is around 368 pg/ml.
Expiration Date	Please enquire.



Mouse
Marapsin/Pancrestin
PicoKine ELISA ...