


Product Datasheet

Human DTK/TYRO3 ELISA Kit (orb381146)

Description	Human DTK/TYRO3 ELISA Kit
Reactivity	Human
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A41-S428
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	96905 MW
Uniprot ID	Q06418
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 1.3-3.5 ng/ml.
Expiration Date	Please enquire.



Human
DTK/TYRO3
PicoKine
ELISA Kit
stand...