


Product Datasheet

Human Fetuin B/FETUB ELISA Kit (orb381147)

Description	Human Fetuin B/FETUB ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C16-P382
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	42055 MW
Uniprot ID	Q9UGM5
Sensitivity	12 pg/ml
Dilution Range	Dilution ratio of 1:2000, concentration in serum and plasma is around 2 ug/ml.
Expiration Date	Please enquire.



Human
Fetuin B
PicoKine
ELISA Kit
standa...