


Product Datasheet

Mouse Persephin ELISA Kit (orb381148)

Description	Mouse Persephin ELISA Kit
Reactivity	Mouse
Range	1.56 pg/ml - 100 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: W22-G156
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	17031 MW
Uniprot ID	070300
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 22 pg/ml.
Expiration Date	Please enquire.



Mouse
Persephin
PicoKine
ELISA Kit
stand...