


Product Datasheet

Human LRIG3 ELISA Kit (orb381152)

Description	Human LRIG3 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D28-T807
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	123434 MW
Uniprot ID	Q6UXM1
Sensitivity	12 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 62.5 pg/ml.
Expiration Date	Please enquire.



Human
LRIG3
PicoKine
ELISA Kit
standard ...