


Product Datasheet

Mouse GDF5/Bmp 14 ELISA Kit (orb381154)

Description	Mouse GDF5/Bmp 14 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A376-R495
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	54895 MW
Uniprot ID	P43027
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 31.2 pg/ml.
Expiration Date	Please enquire.



Mouse GDF5
PicoKine
ELISA Kit
standard c...