


## Product Datasheet

### Human Persephin ELISA Kit (orb381156)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human Persephin ELISA Kit                                                                                                                                                         |
| <b>Reactivity</b>          | Human                                                                                                                                                                             |
| <b>Range</b>               | 93.7 pg/ml - 6,000 pg/ml                                                                                                                                                          |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                             |
| <b>Immunogen</b>           | Expression system for standard: E.coli; Immunogen sequence: W22-G156                                                                                                              |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Sample Types</b>        | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).                                                                                                        |
| <b>MW</b>                  | 16600 MW                                                                                                                                                                          |
| <b>Uniprot ID</b>          | <b>O60542</b>                                                                                                                                                                     |
| <b>Sensitivity</b>         | 10 pg/ml                                                                                                                                                                          |
| <b>Dilution Range</b>      | Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 93.7 pg/ml.                                                                            |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                   |



Human  
Persephin  
PicoKine  
ELISA Kit  
stand...