


Product Datasheet

Mouse TNFRSF25/DR3 ELISA Kit (orb381158)

Description	Mouse TNFRSF25/DR3 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q31-F199
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	44264 MW
Uniprot ID	B1AWN9
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 1000 pg/ml.
Expiration Date	Please enquire.



Mouse
TNFRSF25/DR3
PicoKine ELISA
Kit st...