


Product Datasheet

Human R-Spondin-1 ELISA Kit (orb381159)

Description	Human R-Spondin-1 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S21-A263
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	28959 MW
Uniprot ID	Q2MKA7
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 50-90 pg/ml.
Expiration Date	Please enquire.



Human R-Spondin-1 PicoKine ELISA Kit sta...