


Product Datasheet

Human CD48/Slamf2 ELISA Kit (orb381168)

Description	Human CD48/Slamf2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q27-S220
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27683 MW
Uniprot ID	P09326
Sensitivity	12 pg/ml
Dilution Range	Dilution ratio of 1:20, concentration in serum and plasma is 8-15 ng/ml.
Expiration Date	Please enquire.



Human
CD48
PicoKine
ELISA Kit
standard c...