


Product Datasheet

Human Pan Trypsin ELISA Kit (orb381170)

Description	Human Pan Trypsin ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A16-S247
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	32529 MW
Uniprot ID	P07478
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is around 180 ng/ml.
Expiration Date	Please enquire.



Human Pan
Trypsin
PicoKine
ELISA Kit
sta...