


Product Datasheet

Human Nectin-4/PVRL4 ELISA Kit (orb381173)

Description	Human Nectin-4/PVRL4 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G32-S349
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , cell lysates and urine.
MW	55454 MW
Uniprot ID	Q96NY8
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is 80-500 pg/ml. For plasma, dilution ratio is 1:1 and concentration is around 80 pg/ml.
Expiration Date	Please enquire.



Human
Nectin-4
PicoKine
ELISA Kit
standa...