


Product Datasheet

Mouse Kremen-1 ELISA Kit (orb381174)

Description	Mouse Kremen-1 ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A20-G395
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	51673 MW
Uniprot ID	Q99N43
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 0-80 pg/ml.
Expiration Date	Please enquire.



Mouse
Kremen-1
PicoKine
ELISA Kit
standa...