


Product Datasheet

Mouse LDLR / OLR1 ELISA Kit (orb381176)

Description	Mouse LDLR / OLR1 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A22-R790
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin) and urine.
MW	94947 MW
Uniprot ID	P35951
Sensitivity	12 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is 26.2-62.2 ng/ml.
Expiration Date	Please enquire.



Mouse LDLR
PicoKine
ELISA Kit
standard c...