


Product Datasheet

Human SAA/SAA1 ELISA Kit (orb381177)

Description	Human SAA/SAA1 ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S20-Y122
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	13532 MW
Uniprot ID	PODJ18
Sensitivity	150 pg/ml
Dilution Range	For serum, dilution ratio is 1:20 and concentration is around 140-287 ng/ml. For plasma, dilution ratio is 1:1 and concentration is around 20 ng/ml.
Expiration Date	Please enquire.



Human
SAA/SAA1
PicoKine
ELISA Kit
standa...