


Product Datasheet

Mouse Syndecan-1/SDC1 ELISA Kit (orb381183)

Description	Mouse Syndecan-1/SDC1 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q18-E252
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	32905 MW
Uniprot ID	P18828
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 2 ng/ml.
Expiration Date	Please enquire.



Mouse
Syndecan-
1/SDC1
PicoKine
ELISA Kit...