


Product Datasheet

Mouse Syndecan-3/SDC3 ELISA Kit (orb381185)

Description	Mouse Syndecan-3/SDC3 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A45-V387
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	46002 MW
Uniprot ID	Q64519
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 156 pg/ml.
Expiration Date	Please enquire.



Mouse
Syndecan-
3/SDC3
PicoKine
ELISA Kit...