


Product Datasheet

Rat TNFSF11/RANKL ELISA Kit (orb381188)

Description	Rat TNFSF11/RANKL ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: F141-D318
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	35370 MW
Uniprot ID	Q9ESE2
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 20-200 pg/ml.
Expiration Date	Please enquire.



Rat
TNFSF11/RANKL
PicoKine ELISA
Kit sta...