


Product Datasheet

Rat Sclerostin/SOST ELISA Kit (orb381192)

Description	Rat Sclerostin/SOST ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F29-Y213
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	23974 MW
Uniprot ID	Q99P67
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 400 pg/ml.
Expiration Date	Please enquire.



Rat
Sclerostin/SOST
PicoKine ELISA
Kit s...