

Product Datasheet

Mouse S100A13 ELISA Kit (orb381196)

Description	Mouse S100A13 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-K98
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11158 MW
Uniprot ID	P97352
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is 10-15.9 ng/ml.
Expiration Date	Please enquire.



Mouse
S100A13
PicoKine
ELISA Kit
standar...